



INDUSTRIAL PHARMACY

G D MEMORIAL COLLEGE OF PHARMACY

S. No.	Experiment	Date	Page No.	Marks	Teacher's Signature
1.	To perform <u>preformulation</u> study of given drug - [Paracetamol / aspirin]	08/07/22	01 to 06	9	<i>Am</i>
2.	To prepare Cold cream (20gm)	15/07/22	07 to 08	9	<i>Am</i>
3.	To prepare Vanishing cream (20gm)	22/07/22	09 to 10	9	<i>Am</i>
4.	To prepare Calcium gluconate injection (100ml)	29/07/22	11 to 12	9	
5.	To prepare & submit eye drop & eye ointment	05/08/22	13 to 16	9	
6.	To perform quality control test of marketed tablets as per IP -	09/09/22	17 to 20	08	
7.	To prepare ascorbic acid injection -	16/09/22	21 to 23	09	<i>Am</i>
8.	To preparation & evaluation of tetracycline capsules	30/09/22	24 to 28	09	
9.	To preparation & evaluation of aspirin tablets -	04/10/22	29 to 31	09	
10.	To preparation & evaluation of paracetamol tablets -	01/10/22	32 to 34	08	
11.	To evaluation of glass container as per IP	15/10/22	35 to 37		

G D MEMORIAL COLLEGE OF PHARMACY

S. No.	Experiment	Date	Page No.	Marks	Teacher's Signature
	<p>* Total number of practical = 11</p> <p>* Total marks of practical = 110</p> <p>* Total number of practical in present \Rightarrow</p> <p>* Total obtained marks =</p> $\left\{ \# \text{ Average marks} = \frac{\text{Total obtained marks}}{\text{Total number of practical}} \right\}$ $= 8$				

⇒ Object:- To perform preformulation study of given drug. [Paracetamol / aspirin]-

⇒ Reference: DR. SHEIKH A. ALJAZIR., MD RAGEEB MD USMAN, DR. BIYANJ R. KAILASH "practical book of industrial pharmacy-1 exclusively marketed & distributed by PV books; edition: 2018 & page number: 01 to 09 -

⇒ Requirement:- a). Chemicals: Paracetamol powder & aspirin powder-

b). Apparatus:- Test tubes, pH meter, Funnel, graduated measuring cylinder, stop watch, bulk density apparatus, weighing balance, burette stand, Thiel's tube, thermometer, capillary tube, thread, weighing bottle with stopper, oven, & desiccator.-

⇒ Theory:- → "Preformulation" ←

It can be defined as an investigation of physical & chemical properties of a drug substance only & when combined with ~~expi~~ excipients.

* When a newly synthesized drugs shows a sufficient pharmacologic promise in animal model to warrant evaluation in man.-

⇒ Approximate solubility ⇐

Statement of approximate solubility	Approximate volume of solvent in milliliters per gram of solute
Very soluble	Less than 1
Freely soluble	From 1 to 10
Soluble	From 10 to 30
Sparingly soluble	From 30 to 100
Slightly soluble	From 100 to 1000
Very slightly soluble	From 1000 to 10,000
Insoluble or practically insoluble	More than 10,000

Date. 08/07/22.

S. No. 01

Page No. 02

Procedure :-

- Organoleptic properties:- Check the drug for colour by observing through naked eyes.
 - Take the drug in petridish & observe keenly in day light besides the white & black background.
 - Colour is generally a function of a drug's inherent chemical structure relating to a certain level of unsaturation.
 - Some compound may appear to have colour although structurally saturated.
 - Odor greatly affects the flavor of a preparation or food stuff.
- Identification test:- To a glass test tube added 100mg of drug sample & 10ml of water then added ferric chloride solution. No colour change, remained clear indicate the presence of aspirin (acetylsalicylic acid).
- Solubility:- Solubility is a physical property of a substance in a particular solvent. It is a necessary parameter to check the drug solubility in various solvents before going to actual formulation step.

⇒ Indication of powder flow properties ⇐

Flow- Property	Parameters		
	Carver's index	Hausner's ratio	Angle of repose
Excellent	< 10	1 - 1.11	25 - 30
Good	11 - 15	1.12 - 1.18	31 - 35
Fair	16 - 20	1.19 - 1.25	36 - 40
Passable	21 - 25	1.26 - 1.34	41 - 45
Poor	26 - 31	1.35 - 1.45	46 - 55
Very poor	32 - 37	1.46 - 1.59	56 - 65
Very-very poor	> 38	> 1.60	> 66

Date.. 02/07/2022..

S. No. 01

Page No. 03

- * As solubility of medicament is directly related with the absorption of drug. - II
- * Drug having good solubility dissolved better in the absorption fluid resulting in therapeutic response of it.
- * Examination of solubility is an important prerequisite before actual formulation. -

ex. pH :- The pH of a solution is very important in the field of solubilisation & stability of drugs, as it must be adjusted in formulation for maximum stability.

- * The solution of pH = 7 ⇒ Neutral solution -
- * The solution of pH less than 7 = Acidic solution -
- * The solution of pH more than 7 = Basic solution -
- * The ionic & non-ionic forms of drug depend on the pH of the solution. -
- * For measurement of pH, in a glass beaker prepared aqueous solution of drug & measure the pH of solution.

ex. Bulk density :- It is the weight of powder or granule divided by its volume. -

- * Bulk density is used to check the uniformity in bulk powdered materials, to decide the size of container. -

⇒ Calculation:-

a) Bulk density - Mass = 17 gm -
 Volume = 38 ml -
 Bulk density (gm/ml) = $\frac{\text{Weight of sample (W)}}{\text{Bulk volume (V)}}$
 $= \frac{17}{38} = \underline{\underline{0.44 \text{ gm/ml}}}$

b) Tapped density - Mass = 17 gm -
 Volume = 29 ml -
 Tapped density (gm/ml) = $\frac{\text{Weight of sample (W)}}{\text{Tapped volume (V)}}$
 $= \frac{17}{29} = \underline{\underline{0.58 \text{ gm/ml}}}$

c) Carri's index -

Carri's index (%) = $\frac{\text{Tapped density} - \text{Bulk density}}{\text{Tapped density}} \times 100$
 $= \frac{0.58 - 0.44}{0.58} \times 100$
 $= \frac{0.14}{0.58} \times 100, = \underline{\underline{24.13\%}}$

Date..08/07/22.

S. No....01.....

Page No.....04.....

- Weigh accurately 17 gm of drug powder (W). -
- Place it in dried graduated measuring cylinder & note the volume occupied, apparatus for 3 tapping. (V)
- Calculate the bulk density by using formula -

f) Tapped density:- Weigh accurately 17 gm of drug powder (W). -
 - Place it in dried graduated measuring cylinder & note the ^{tapped} volume, apparatus for 100 tapping
 - Calculate the tapped density by using formula -

g) Carri's index:- The Carri's index is an indication of the compressibility of a powder. -
 - This index is frequently used in pharmaceuticals as an indication of the flowability of a powder. -

h) Hausner's ratio:- It is a number that is correlated to the flow-ability of a powder material. -
 - It is used in wide variety of industries as an indication of flow-ability of a powder. -

i) Angle of repose:- The flow properties of powders are critical for an efficient tableting operation.

d). Hausner's ratio -

$$\text{Hausner's ratio} = \frac{\text{Tapped density}}{\text{Bulk density}}$$

$$= \frac{0.58}{0.44} = \underline{\underline{1.31}}$$

e).

Angle of repose -

$$\tan \theta = \frac{h}{r}$$

where:

θ = angle of repose -

h = height of pile (cm) -

r = radius of circle (cm) -

$$\theta = \tan^{-1}\left(\frac{h}{r}\right)$$

$$\left\{ \begin{array}{l} h = 4 \text{ cm} \\ r = 3.7 \end{array} \right\}$$

$$\theta = \tan^{-1}\left(\frac{4}{3.7}\right) = \tan^{-1}(1.08)$$

$$\underline{\underline{\theta = 47.2^\circ}}$$

Date. 02/07/2022.

S. No. 01

Page No. 05

* A good flow of the powder to be compressed is necessary to assure efficient mixing & acceptable weight uniformity for the compressed tablets -

- Take a dry & clean funnel with a round stem of 20-30mm diameter with flat tip -
- Attach the funnel to the burette stand & on a clean, dry platform place a graph paper sheet below the funnel -
- Adjust the distance of 4cm b/w the lower tip of funnel & sheet -
- Pour the sample gently in funnel so that a heap of powder forms & touches the lower tip of the funnel -
- Draw the circle around the heap by using pencil -
- Measure radius & height of the cone & put value in equation to get the value of angle of repose -

j). Melting point :- It is the temperature at which a solid becomes a liquid at normal atmospheric pressure -

- It is an identity for pure crystalline organic compounds -
- Take a capillary tube & at one end seal it by applying flame for a very short time -
- Added a small quantity of drug powder (3-5mg) from another side.

⇒ Observation table for preformulation study on aspirin drug

Test	Observation
Colour -	White crystalline powder
Odour -	Odorless
Taste -	Bitter taste
Identification taste -	No colour change
Solubility with water -	Slightly soluble (0.33 gm/100ml)
pH -	3.5
Bulk density -	0.44 gm/ml
Tapped density -	0.58 gm/ml
Hausner index -	24.13%
Hausner ratio -	1.31
Angle of repose -	47.2°
Melting point -	138°C

Date.. 08/07/22

S. No. 01

Page No. 06

- Tie the capillary tube containing sample with thermo-meter near its bulb with thread.
- Added the liquid paraffin in Thiele's tube & immerse the thermo-meter with capillary tube into liquid paraffin.
- Care should be taken so that thread should not come in contact with liquid paraffin & rate of heating should not be more than 3°C/minute.
- Apply heat to the Thiele's tube with the help of burner & carefully record the temp. at which solid drug compound converts to liquid.
- Repeat the procedure twice & note down the readings carefully.

⇒ Result:-

Preformulation studies of aspirin drug was observed completed.

⇒ Object:- To prepare of cold cream (20gm) -

⇒ Reference:- "Gand R.S. & Gupta G.D." Pharmaceuticals practical book published by Satish kumar Jain, 1st edition : 2002 ; Page number : 32 -

⇒ Requirement - (A) Chemicals:- Cetyl esters wax, white wax, mineral oil,

Sodium borate etc.-

(B) Apparatus:- Beaker, measuring cylinder, glass stirring rod, water bath & weigh balance etc.-

⇒ Theory:- Cold cream is a water-in-oil emulsions.-
- Creams are semisolid dosage forms.-

* Application of cold cream avoids aggravation of skin problems during the cold season.-

* Cold cream is primarily used to prevent excessive drying of skin.-

* It has emollient action & prevents dehydration of skin.-

* The key chemical constituent of cold cream is bees wax.-

* Cold cream is mainly used for skin treatment, due to its moisturizing properties.-

⇒ Formula - (20gm) -

Ingredients	Quantity given (gm)	Quantity taken (gm)
• Cetyl esters wax -	12.5 gm	2.5 gm
• White wax -	12.0 gm	2.4 gm
• Mineral oil -	56.0 gm	11.2 gm
• Sodium borate -	0.5 gm	0.1 gm
• Purified water	19.8 gm	3.8 gm
	Total = 100 gm	Total = 20 gm

⇒ Calculation -

$$\left\{ \begin{array}{l} \text{Taken quantity (x) gm} = \frac{\text{Given quantity}}{\text{Total amount}} \times \text{Preparation amount} \end{array} \right\}$$

eg. - Cetyl esters wax = $\frac{12.5}{100} \times 20 \Rightarrow 2.5 \text{ gm}$

Date... 15/07/22...

S. No... 02

Page No... 08

⇒ Procedure - Reduce size of cetyl esters wax & white wax in small pieces. -

- Melt them on a steam bath with the mineral oil & continue heating until temp. of the mixture reached 70°C (A). -

- Dissolve the sodium borate in purified water & heat at 70°C in separate flask -

- Gradually added aqueous solution (B) to the mixture of oily phase (A) with continues & rapid stirring until it has congealed. -

⇒ Category - Emollient & cleansing cream. -

⇒ Storage - Preserve in a well-closed container. -

⇒ Instruction - For external use only. -

- Not for ophthalmic, oral, or intra-vaginal use. -

⇒ Result - Cold cream was prepared & submitted. -

⇒ Object:- To preparation of vanishing cream (20gm) -

⇒ Reference:- "Gand R.S. & Gupta G.D." pharmaceuticals practical book published by Satish Kumar Jain, 1st edition: 2002; page number: 33 -

(A). Chemicals -

⇒ Requirements:- Stearic acid, Stearyl alcohol, Cetyl alcohol, methyl paraben, propyl-paraben, potassium hydroxide & purified water etc -

(B). Apparatus -

- Beaker, measuring cylinder, glass stirring rod, water bath & weigh balances etc -

⇒ Theory:-

A cream is a preparation usually for application to the skin -

- Creams are semi-solid emulsions of oil & water -
- Vanishing creams get their name from the fact that they seem to disappear when spread on the skin -
- It is also known as foundation cream which are applied to skin to provide a smooth emollient base before the application of face powder & other face make up -

* Formulation - [Formula] - (20gm) -

Ingredients	Quantity given (gm)	Quantity taken (gm)
• Stearic acid -	13 gm.	2.6 gm -
• Stearyl alcohol -	1.0 gm -	0.2 gm -
• Cetyl alcohol -	1.0 gm -	0.2 gm -
• KOH (potassium hydroxide) -	0.90 gm.	0.18 gm -
• Methyl paraben -	0.10 gm -	0.02 gm
• Propyl paraben -	0.05 gm -	0.01 gm
• Purified water - (q.s.)	83.95 gm.	16.79 gm

* Formula $\Rightarrow \frac{\text{Given quantity}}{\text{Total amount}} \times \text{Prepare quantity} -$

Date. 12/07/22...

S. No. 03.

Page No. 10

- Chemically they are oil in water emulsions consisting of stearic acid, an alkali, a polyol & water.
- One important constituent of the cream is "hydroquinone".

\Rightarrow Procedure -

- melt stearic acid, stearyl alcohol & cetyl alcohol at 75°C in a conical flask. (A).
- Dissolve potassium hydroxide in purified water at 75°C in another conical flask & added preservative. (B).
- Added (A) into (B) at 75°C slowly with continuous stirring.
- slowly stir the mixture until a smooth cream was formed at room temperature.

\Rightarrow Category :- Vanishing cream -

\Rightarrow Storage :- store in well closed container, in cool place. Do not freeze.

\Rightarrow Instruction - for external use only. Not for Ophthalmic, Oral, or intravaginal use.

\Rightarrow Result :- Vanishing cream was prepared & Submitted.

* Object:- To prepare Calcium gluconate injection [100ml]-

* Reference- "Gand R.S. & Gupta G.D." pharmaceutical book published by Satish Kumar Jain, 1st edition: 2002; page number: 153 to 154.-

* Requirement - (A) Chemicals-

• Calcium gluconate, Calcium-D-saccharate & water for injection etc.-

(B) Apparatus- Beaker, measuring cylinder, syringe, needle, ampoules, Whatman filter paper / membrane filter, autoclave.-

* Theory:- Calcium is the most abundant cation & the 5th most common inorganic element in the human body.-

* Calcium is essential for the maintenance of the nervous muscular & skeletal systems, & for cell membrane & capillary permeability.-

* Its role in bone structure & muscle contraction is well known, but Ca^{2+} is also important for blood

⇒ Formula -

Ingredient	Given quantity	Taken quantity
1. Calcium gluconate -	5.0 gm	5.0 gm
2. Calcium-D-saccharate-	1.75 gm	1.75 gm
3. Water for injection	Q.S. to 100ml	100ml (Q.S.)

⇒ Taken quantity =

[Formula] ⇒ $\frac{\text{Given quantity}}{\text{Total amount}} \times \text{Prepare quantity}$

eg. Calcium gluconate ⇒ $\frac{5.0}{100} \times 100 \Rightarrow \underline{5.0 \text{ gm}}$

• Calcium-D-saccharate ⇒ $\frac{1.75}{100} \times 100 \Rightarrow \underline{1.75 \text{ gm}}$

• Water for injection ⇒ $\frac{100 (Q.S.)}{100} \times 100 \Rightarrow \underline{100 \text{ ml (Q.S.)}}$

Date..29/05/2022

S. No. 04

Page No. 12

Cogulation, nerve conduction & electrical conduction in the myocardium.-

* Calcium gluconate is calcium salt of gluconic acid, an oxidation product of glucose.-

* Procedure - All glassware were sterilized by the autoclave at 121°C for 15 minutes.-

- Weigh accurately calcium gluconate & Calcium D-saccharate.-
- Dissolve calcium gluconate in water for injection in a beaker with application of heat then added Calcium D-saccharate to it.-
- Adjust the pH in between 6-8 with 10% NaOH solution.-
- Keep the solution for cooling.-
- After cooling filter the ~~sol~~ solution through Whatman filter paper / 0.45 µm membrane filter to remove any particulate matter.-
- Fill the prepared injection in conical flask with the help of measuring cylinder.-
- Seal the conical flask by aluminium foil paper.-
- Sterilized the conical flask by autoclaving at 121°C for 30 minutes.-

* Result :- Calcium gluconate injection was prepared & submitted successfully.-

* Object - To prepare & submit eye drop & eye ointment -

Reference - "Gand R.S. & Gupta G.D." pharmaceuticals practical book published by Satish Kumar Jain, 1st edition: 2002, page number: 130, 131 & 142 -

Requirement -

(A) 'Chemicals' -

i) Solution for eye drops -

• Methyl hydroxybenzoate,
Propyl hydroxybenzoate, Purified water -

ii) Zinc sulphate eye drop -

• Zinc sulphate, Sodium chloride & solution for eye drops q.s. to -

iii) Eye ointment -

• Woolfat, Yellow soft paraffin & Liquid paraffin q.s. to -

(B) Glassware -

• Beaker, glass stirring rod, measuring cylinder, spatula, heating mantle, autoclave & aseptic rooms -

Formula -

① Solution for eye drops - (50ml)

"Ingredient"	"Given quantity"	"Taken quantity"
① Methyl hydroxy benzoate	22.0 mg	11.0 mg
② Propyl hydroxy benzoate	11.4 mg	5.7 mg
③ Purified water, sufficient to produce -	100.0 ml (q.s.)	50 ml (q.s.)

Taken quantity = $\frac{\text{Given quantity}}{\text{Total amount}} \times \text{Prepare quantity}$

eg. Methyl hydroxy benzoate = $\frac{22}{100} \times 50 = 11.0 \text{ mg}$ -

Propyl hydroxy benzoate = $\frac{11.4}{100} \times 50 = 5.7 \text{ mg}$ -

Purified water (q.s.) = $\frac{100}{100} \times 50 \Rightarrow 50 \text{ ml (q.s.)}$ -

Date..05/08/22

S. No. 05

Page No. 14

Theory:- i. Solution for eye drops - (50ml) -

- Solution for eye drops must be freshly prepared aseptically & dispensed in previously sterilized containers.

ii. Zinc sulphate eye drop - (10ml) -

- Ophthalmic preparations are sterile products, free from foreign particles.
- Aqueous eye drops contain suitable anti-microbial preservatives.

* Eye drops are sterile aqueous / oily solutions / suspension.

① Examples of eye drops - Zinc sulphate eye drop, Framycetin eye drop, Pilocarpine eye drop, Silver nitrate eye drop.

iii. Eye ointment - (20gm)

- Eye ointments must be free from large particles & should be prepared under aseptic condition.

- Ophthalmic ointments are help to improve ocular bioavailability & sustain the drug release.

① Example of eye ointment - Chloramphenicol ointment, tetracycline ointment & Hydro-cortisone ointment.

Formula -

① Zinc sulphate eye drop - (10ml)

"Ingredient"	"Given quantity"	Taken quantity
① Zinc sulphate -	12.5 mg -	25.0mg
② Sodium chloride -	40.0 mg -	80.0mg
③ solution for eye drops q.s. to	5.0ml (q.s.)	10.0ml (q.s.) -

Taken quantity = $\frac{\text{Given quantity}}{\text{Total amount}} \times \text{Prepare quantity}$

eg. zinc sulphate = $\frac{12.5}{5} \times 10 \Rightarrow 25.0\text{mg}$ -

• Sodium chloride = $\frac{40}{5} \times 10 \Rightarrow 80.0\text{mg}$ -

• solution for eye drops (q.s.) = $\frac{5}{5} \times 10 \Rightarrow 10\text{ml (q.s.)}$ -

Date. 05/08/2022

S. No. 05

Page No. 15

Procedure - All glass-ware were sterilized by the auto-clave at 121°C for 15 minutes. -

↓

i. Solution for eye drops - (50ml) -

- Dissolve methyl hydroxy-benzoate & propyl hydroxy-benzoate in boiling water under aseptic condition. -
- Added freshly boiled & cooled purified water to produce the required volume. -

ii. Zinc sulphate eye drop - (10ml) -

- Dissolve weighed amount of zinc sulphate in purified water with aseptic condition. -
- Dissolve sodium chloride in small amount of purified water -
- Added this solution & mix it required filter it aseptically condition. -

iii. Eye ointment - (20gm) -

- melt together weighed amount of wool fat & yellow soft paraffin in a container. -
- Added liquid paraffin to make up weight 10gm. -

Formula-

○ "Ointment eye" - (20 gm)

○ "Ingredient"	○ "Given quantity"	Taken quantity
○ Zinc sulphate -		
○ Wool fat	1.0 gm	2.0 gm
○ Yellow soft paraffin -	8.0 gm	16.0 gm
○ Liquid paraffin, sufficient to produce -	q.s. 10 gm -	20 gm (q.s.)

$\text{Taken quantity} = \frac{\text{Given quantity}}{\text{Total amount}} \times \text{Prepare quantity}$

eg. - Wool fat = $\frac{1.0}{10} \times 20 = 2.0 \text{ gm}$

• Yellow soft paraffin = $\frac{8.0}{10} \times 20 = 16.0 \text{ gm}$

• Liquid paraffin (q.s.) = $\frac{10}{10} \times 20 = 20 \text{ gm (q.s.)}$

Date.. 05/08/22..

S. No. 05

Page No. 16

- Filter the hot mixture through coarse filter paper placed in a heated funnel. -
- Filterate was sterilized by dry heat at 150°C for sufficient time to ensure that the whole was maintained at this temperature for one hour. -
- Allow to cool at room temperature without opening the container. -

Results-

- Eye drop & eye ointment were prepared & submitted. -

Object:-

To perform quality control test of marketed tablets as per IP-

Reference:-

Indian pharmacopoeia-2010 published by the indian pharmacopoeia commission Ghaziabad (U.P.), 6th edition effective from 1st september: 2010, page number: 187-193 (Volume-1st)
& 751 - 754 (Volume-2nd)



Requirement:-(A) Chemicals-

• Tablets, 0.1N Hydrochloric acid (HCl), phosphate buffer at pH 6.8 [Potassium dihydrogen phosphate & disodium hydrogen phosphate] etc-

(Basket-rack assembly)

(B) Apparatus-

• Friability apparatus, Disintegration apparatus, volumetric flask, Beaker, measuring cylinder, weigh balances & Whatman filter paper, Monsanto hardness tester / pfizer hardness tester & Vernier caliper etc-



Date 09/09/2022

S. No. 06

Page No. 18

Theory :- # Tablets

- Tablet is a solid unit dosage form.
- The solid unit dosage form of medication/medicament with suitable excipients.
- * Tablets are intended for oral administration.
- Tablets are prepared by the "compression method" are called compressed tablets.
- * Tablets \Rightarrow Circular in shape, Biconvex/flat in shape.

Quality Control test # (Quality) -

* marketed tablet test (Finished product test) -

i). Disintegration test / disintegrating test -

• This test is

not used for "chewable & sustain release" tablets.

- Uncoated tablets \Rightarrow water, UMT - 15 minutes.
 - Film coated tablets \Rightarrow water, UMT - 30 minutes.
 - Coated tablets \Rightarrow water, UMT - 60 minutes.
 - Enteric coated tablets \Rightarrow In 0.1N HCl solution operate "without the discs" for 2 hrs \downarrow
- No tablets shows signs of cracks / no change in tablets.
 - \downarrow
 - At pH 6.8 phosphate buffer solution operate "with discs"
 - \downarrow
 - UMT - 30 minutes.

Description

- 7x4 mm, Round shape, biconvex
 - Yellow to light brown colour
- ⇒ Enteric coated tablets.

* Average weight ⇒ 0.169 gm - (169 mg)

* Diameter = 7 mm.

* Thickness = 4 mm.

* Hardness = 15 kg/cm².

* Friability test = NA

* Disintegration time = In 0.1M HCl solution for 2 hrs
Volume 900ml & temp. 37°C ± 2°C

↓
After 2 hrs no change in tab.
↓

At pH 6.8 phosphate buffer solution, Volume 900ml & 37°C ± 2°C

↓
NMT 60 minutes

Date 09/09/2022

S. No. 06

Page No. 19

v. Dispersible tablets ⇒ 24°C to 26°C, NMT - 3 minutes -

* Uniformity of dispersible tab. ⇒ 2 tablets in 100ml water ⇒ eff. for
2 & 1/2 hrs

vi. Effervescent tablets ⇒ 20°C to 30°C, 1 tablet in 250ml water, NMT - 5 minutes

vii. Soluble tablets (Coated/film coated) ⇒ 15°C to 25°C, NMT - 3 minutes -

ii. Dissolution test -

• The test is studied for a tablet & capsule -

* Apparatus as per IP ⇒ Paddle (1st) & Basket (2nd) - 7.5 rpm, 900ml

* Apparatus as per USP ⇒ Basket (1st) & Paddle (2nd) & 37°C -

iii. Content of active ingredient (Assay) - ± 10%

• 10 tablets cross ⇒ 1 tablet blender ⇒ Dilution & absorption -

iv. Content uniformity test -

• Each individual tablet less than

10mg drugs or < 10% drugs -

v. Hardness - used apparatus - * Monsanto hardness tester -

* Pfizer & Storkab tester -

Procedure -

• Put one tablet into each tube, suspend the assembly in the beaker containing 0.1M HCl & operate without the discs for 2 hrs (hours) -

Observation table

* Diclofenac sodium tablets IP 50mg *

* B.No.: C20RDA048 -

S. No.	Tablets	Weight (gm)	Diameter (cm)	Thickness (cm)	Hardness
1.	Tablet-I	0.176 gm	0.7cm	0.4cm	15 kg/cm ²
2.	Tablet-II	0.174 gm	0.7cm	0.4cm	
3.	Tablet-III	0.171 gm	0.7cm	0.4cm	
4.	Tablet-IV	0.163 gm	0.7cm	0.4cm	
5.	Tablet-V	0.168 gm	0.7cm	0.4cm	
6.	Tablet-VI	0.166 gm	0.7cm	0.4cm	

Date: 09/09/2022

S. No. 06

Page No. 20

- Unless otherwise stated in the individual monograph.
- Remove the assembly from the liquid.
- No tablet shows signs of cracks that would allow the escape of the contents on disintegration, apart from fragments of coating.
- * Replace the liquid in the beaker with mixed phosphate buffer pH 6.8, add a disc to each tube & operate the apparatus for further UMT 60 minutes.

Result

• Quality control test of marketed tablet (Diclofenac sodium tablets IP 50mg) as per IP was observed.

Object:-

To prepare ascorbic acid injection (20ml) -

Reference:-

• "DR. Kailash R. Blyani, MD. Razeed MD. Usman, DR. Aijaza A. Shetkh" practical book of industrial pharmacy-1st published by Pee Vee, Edition: 2018, page number: 59-62- & Gaud & Gupta-Page number-151-

Requirement:-(A) Chemical-

• Ascorbic acid, Disodium EDTA, Sodium hydroxide, Sodium bicarbonate & water for injection q.s. - , propylene glycol-

(B) Apparatus-

• Beaker, measuring cylinder, syringe, needle, amber colour ampoules, membrane filter, pH meter & weigh balance etc-

Theory-# Ascorbic acid

- Ascorbic acid (Vitamin-C) is a water-soluble vitamin.-
- It occurs as a white or slightly yellow crystal or powder with a light acidic taste.-

formula -

• Ascorbic acid injection - (20ml) -

Ingredients.	Given quantity.	Taken quantity.
• Ascorbic acid.	1.03 gm.	0.41 gm.
• Sod. bicarbonate.	2.47 gm.	0.98 gm.
• Propylene glycol.	2.50 gm.	1.00 gm.
• Sod. hydrosulphite.	0.25 gm.	0.1 gm.
• Disod. edetate.	0.05 gm.	0.02 gm.
• Water for injection (q.s.)	50.0 ml (q.s.)	20 ml (q.s.)

$$* \text{Taken quantity} = \frac{\text{Given quantity}}{\text{Total amount}} \times \text{Prepare quantity/amount}.$$

$$* \text{Ascorbic acid} = \frac{1.03}{50} \times 20 = 0.41 \text{ gm.}$$

Date... 16/09/22

S. No... 07

Page No... 22

- It is an anti-scorbutic product.
- The chemical name of ascorbic acid is L-ascorbic acid.
- * Ascorbic acid injection is a clear, colorless to slightly yellow sterile solution of ascorbic acid in water for injection for i.v., i.m. & subcutaneous used & contains no preservatives.

Procedure -

- The formulation must be carried out in clean area.
- Weigh the ascorbic acid accurately.
- Dissolve ascorbic acid slowly under constant stirring & nitrogen bubbling.
- Added sod. bicarbonate & slow stirring.
- Separately prepare solution of sodium hydrosulphite in small proportion of water for injection & added in the solution.
- Added solution of disodium edetate & mix uniformly.
- Added mix propylene glycol with slow stirring.
- pH is adjusted to 5.5 to 6.5 using sod. hydroxide solution or ascorbic acid solution.
- make up the volume & mix with water for injection.

• Sodium bicarbonate = $\frac{2.97}{50} \times 20 = 0.988 \text{ gm} -$

• Propylene glycol = $\frac{2.50}{50} \times 20 = 1.00 \text{ gm} -$

• Sodium hydrosulphite = $\frac{0.25}{50} \times 20 = 0.1 \text{ gm} -$

• Disodium edetate = $\frac{0.05}{50} \times 20 = 0.02 \text{ gm} -$

• Water for injection (q.s.) = $\frac{50}{50} \times 20 = 20 \text{ ml (q.s.)} -$

Date... 16/09/2022

S. No... 07

Page No... 23

• Sterilized by passing through sterilized 0.22 μ membrane filter into sterilized filling vessel. -

Result -

• Ascorbic acid injection ampoules were prepared & submitted. -

Object :-

- Preparation & evaluation of capsules -
(Tetracycline capsules = 250mg) -

Reference -

- DR. Shekh Aijaz, DR. Kallash "Industrial pharmacy - 1st practical book published by Pee Dee, Edition: 2018 & page number - 43 to 52."


Requirement -

- Chemical - Tetracycline hydrochloride, dried starch & dried talcum
- Apparatus -
• Mortar & pestle, empty capsules, sieve, hand operated capsule filling machine, volumetric flask, pipette, beaker, stop watch, measuring cylinder, Whatman filter paper & Disintegration test apparatus.

Theory -"CAPSULES"

- Capsules are solid dosage form in which drugs is enclosed in hard gelatin capsule shell -
- in soft soluble shell of gelatin -
- * Various shapes & Capacities -

Tetra-cycline



- Colour - Yellow -
- Odour - Odourless -
- * Crystalline powder -
- Soluble - water, acid & alkali -
- Insoluble - Chloroform & ether -

* It is stable in air -

• But exposure to strong sunlight causes it to darken -

* These are readily absorbed & bound to plasma protein in varying degrees -

Uses - It is a prescription antibiotic -

- Treat a wide range of infections -
- eg - Skin infection -
- Respiratory tract infection -
- Genital & urinary infection -

S. No. 08

Date 30/09/2022

Page No. 25

Preparation of tetracycline Capsules -

Procedure

- Weigh the required quantity of drug & excipients -
- Moisture content should be less than 1.5% -
- starch should be dried & sieved through # 100, moisture less than 1.5% -
- Talcum should be dried & sieved through # 100, moisture less than 1% -
- mix all the ingredients uniformly using mortar & pestle -
- Empty capsule shells of number "0" is selected & for filling the content in capsule shell hand operated capsule filling machine is used -
- ↓
- This machine consist of a bed, a loading tray, a powder tray, a pin plate, a sealing plate having a rubber top, a lever, a cam handle -
- ↓
- The empty capsules are filled in the loading tray & it is placed over the bed -
- ↓
- The cam handle is operated to separate the Capsule caps from their bodies -
- ↓

formula-

Ingredient	Given quantity	Taken quantity
• Tetracycline Hcl -	250 mg -	208.33 mg -
• Dried starch -	25 mg -	20.83 mg -
• Dried talcum -	25 mg -	20.83 mg -

$$\# \cdot \text{Taken quantity} = \frac{\text{Given quantity}}{\text{Total quantity}} \times \text{Prepare quantity}$$

$$\cdot \text{Tetracycline Hcl} = \frac{250}{300} \times 250 = 208.33 \text{ mg}$$

$$\cdot \text{Dried starch} = \frac{25}{300} \times 250 = 20.83 \text{ mg}$$

$$\cdot \text{Dried talcum} = \frac{25}{300} \times 250 = 20.83 \text{ mg}$$

Date. 30/09/22.

S. No. 08

Page No. 26

- The powder tray is placed in a proper position & filled with an accurate quantity of powder with scraper.
⇓
- The excess of the powder is collected on the platform of the powder tray -
⇓
- The pin plate is lowered & the filled powder is pressed by moving the pin downwards.
⇓
- After pressing the pin plate was raised & the remaining powder was filled into the bodies of the capsules
⇓
- The powdered tray was removed after its complete filling -
⇓
- The cap holding tray was again placed in position -
⇓
- The plate with the rubber top was lowered & the lever was operated to lock the caps & bodies -
⇓
- The loading tray was then removed & filled capsules are collected. -
⇓
- Capsules were store in air tight container -

Evaluation of tetracycline Capsules-i) Weight variation/uniformity of weight test-

- 20 capsules are used individually weigh & calculate avg.
- Average weight of capsule 0% deviation/variation ^{net-}

Contents

$$\cdot < 300\text{mg} = 10\%$$

$$\cdot \geq 300\text{mg} = 7.5\%$$

ii) Disintegration test - (H4C) - • Water (medium) -

- Operate the apparatus for 30 minutes.

iii) Content of active ingredients - (Assay) - $\pm 10\%$ -

- This range is based on the requirement that 20 capsules are used in the assay.

* NLT 5 Capsules are used.

iv) Uniformity of Content -

- This test is applicable to capsules that active ingredient contain less than 10mg/ 10% w/w -

* For Capsules -

• Containing more than one active ingredient

- Test for each active ingredient that corresponds to the afore-mentioned conditions -

* 10 Capsules taken at random -

- NMT 1- individual values outside the 85% - 115% -
- Not any individual values outside the 75% - 125% -



- If 2/3 individual value are outside 85% - 115% -



- 20- capsules another taken (Total = 30 Capsules (10+20)) -
- NMT 3- individual values outside the 85% - 115% -
- Not any individual values outside 75% - 125% -

Result :-

• Tetracycline capsules were prepared & evaluated -

Object-

To preparation & evaluation of aspirin tab. -
(75mg - 20 tabs)

Reference-

"DR. Shekh Aijaz, DR. Kailash" industrial pharmacy - 1st practical book published by Pee Dee.
Edition: 2018 & page number: 23 - 32 -

Requirement -

(A) Chemical -

Aspirin, micro crystalline cellulose (mcc - Avicel PH-102), Croscarmellose sodium, Talc & magnesium stearate -

(B) Apparatus -

Single punch tablet machine, weight balance, mortar & pestle, spatule & Sieves #40, #44, #60 & #100 -

Theory -

Aspirin or acetylsalicylic acid is a derivative of salicylic acid -

Aspirin is a non-steroidal anti-inflammatory drugs (NSAID) -

Formula - [1 tab. = 150mg, 20 tab x 150 = 3000mg]

Ingredient	Given quantity	Taken quantity
Aspirin-	75mg-	1.5gm-
microcrystalline cellulose (MCC)-	52.5mg-	1.05gm-
Croscarmellose sodium-	7.5mg-	0.15gm-
Talc-	7.5mg-	0.15gm-
magnesium stearate-	7.5mg-	0.15gm-

• Taken quantity = $\frac{\text{Given quantity}}{\text{Total amount}} \times \text{Preparation amount}$

i. Aspirin = $\frac{75}{150} \times 3000 = 1500\text{mg} / 1.5\text{gm}$

ii. MCC $\Rightarrow \frac{52.5}{150} \times 3000 = 1050\text{mg} / 1.05\text{gm}$

iii. Croscarmellose sodium = $\frac{7.5}{150} \times 3000 = 150\text{mg} / 0.15\text{gm}$

iv. Talc = $\frac{7.5}{150} \times 3000 = 150\text{mg} / 0.15\text{gm}$

v. Magnesium stearate = $\frac{7.5}{150} \times 3000 = 150\text{mg} / 0.15\text{gm}$

Date.. 01/10/2022

S. No... 09

Page No... 30

Procedure - [By direct compression method] -

- Weigh the required quantity of aspirin & other excipients.
- Pass the aspirin & other excipients through #40/44 & mixed for 15 minutes.
- mixed thoroughly all the excipients with granular aspirin in mortar with the help of spatula in order to get uniform distribution of ingredients.
- Pass the talc & magnesium stearate through #60 - mixed uniform distribution of ingredients for 2 minutes.
- The granules are ready for compression & can be compressed into tablets.

Evaluation of aspirin tablets -

i. Weight variation - 20 tablets used for test -

Weight

% Variation

○ $\leq 80\text{mg}$ -



10%

○ $>80\text{mg} - <250\text{mg}$ -

7.5%

○ $\geq 250\text{mg}$ -

5%



Date. 01/10/22.

S. No. 09

Page No. 31


- ② Friability test - $\text{Weight} \geq 6.5 \text{ gm} / 7 \text{ gm}$ -
 - 25 rpm & 4 minutes \Rightarrow 100 rotation -
 • $\text{Weight} \Rightarrow$ weight loss NMT 1% -

- ③ Dissintegration time -
 • Aspirin uncoated tablets -
 • Uncoated tablets -
 water, NMT - 15 minutes -

- ④ Hardness -
 • The hardness of a tablet is an indication of its strength.
 - Depending upon the type & conc. of the binding agent the hardness of tablet varies -
 * The force is measure in kg/cm^2 -

- # Result -
 • Aspirin tablets was prepared & evaluated -

Paracetamol # [Acetaminophen-]



* White or almost white -

* Crystalline powder -

* Standard (IP) -

• PCM tablets contain not less than 95% & not more than 105% of the stated amount of PCM ($C_8H_9NO_2$) -

Date. 01/10/22.

S. No. 10

Page No. 32

Object -

• To preparation & evaluation of paracetamol tablets (500mg) -

Reference -

• "SR. Shekh Aijaz, SR. Kallash" industrial pharmacy - 1st practical book published by Pee Vee. Edition: 2018 & page number: 12 to 22 -

Requirement -

(A) Chemicals -

• Paracetamol, microcrystalline cellulose (mcc), Croscarmellose sodium, Poly-vinyl-pyrrolidone (PVP), Talc, magnesium stearate & purified water -

(B) Apparatus -

• Weigh balance, single punch tablet machine, mortar & pestle, measuring cylinder, beaker, Sieves #40/44 & #60 - & dryer -

Theory -

• Paracetamol was discovered in 1877 -
• PCM may inhibit the nitric oxide (NO) pathway mediated by a variety of neurotransmitter receptors -

Formula - [1 tab - 800 mg, 20 tabs = 16,000 mg] -

Ingredient	Given quantity	Taken quantity
Paracetamol -	500 mg -	10 gm -
microcrystalline cellulose -	180 mg -	3.6 gm -
croscarmellose sodium -	40 mg -	0.8 gm -
Polyvinyl pyrrolidone (K-30)	40 mg -	0.8 gm -
Talc	20 mg -	0.4 gm -
magnesium stearate -	20 mg -	0.4 gm -
Purified water	Q.S. (10ml)	Q.S. (10ml) -

• Taken quantity = $\frac{\text{Given quantity}}{\text{Total amount}} \times \text{Prepare amount}$ -

i. PCM = $\frac{500}{800} \times 16,000 = 10,000 \text{ mg} \mid 10 \text{ gm}$ -

ii. MCC = $\frac{180}{800} \times 16,000 = 3,600 \text{ mg} \mid 3.6 \text{ gm}$ -

iii. Croscarmellose sod. = $\frac{40}{800} \times 16,000 = 800 \text{ mg} \mid 0.8 \text{ gm}$ -

iv. PVP = $\frac{40}{800} \times 16,000 = 800 \text{ mg} \mid 0.8 \text{ gm}$ -

v. Talc = $\frac{20}{800} \times 16,000 = 400 \text{ mg} \mid 0.4 \text{ gm}$ -

vi. mag. stearate = $\frac{20}{800} \times 16,000 = 400 \text{ mg} \mid 0.4 \text{ gm}$ -

Date.. 01/10/22..

S. No. 10

Page No. 33

Procedure - [By wet granulation] -

- 1st • Weigh the required quantity of Paracetamol & other all excipient -
- 2nd. Paracetamol, microcrystalline cellulose & croscarmellose sodium were passed through sieve # 40/44 & mixed for 15 minutes -
- 3rd. Polyvinyl pyrrolidone (K-30) was dissolved in 10 ml purified water -
- 4th. Step-2nd materials were granulated with step-3rd binder solution with "Continuous kneading" -
- 5th. Step-4th materials were passed through sieve # 10 -
- 6th. Step-5th materials was dry at 60°C for 20 minutes in dryer -
- 7th. Step-6th materials was passed through sieve # 40/44 -
- 8th. Talc & magnesium stearate was passed through sieve # 60 & mixed with step-7th material for 5 minutes -
- 9th. The granules are ready for compression & can be compressed into tablets -

Weight Evaluation of Paracetamol tab -

① Weight variation - 20 tab. used for test -

Weight	% Variation
- $\leq 80\text{mg}$	10%
- $> 80\text{mg} - < 250\text{mg}$	7.5%
- $\geq 250\text{mg}$	5%

② Friability test - Its integration time -

- ~~PCM~~ uncoated tablets -
- uncoated tab.
- Water, NMT - 15 minutes -

③ Friability test -

- Weight $\geq 6.5\text{gm} / 7\text{gm}$
- 25 rpm \times 4 minutes \Rightarrow 100 rotation -
- Weight \Rightarrow weight loss NMT 1% -

④ Hardness - The hardness of a tablets is an indication of its strength. -

- Depending upon the type & Conc. of the binding agent the hardness of tablet varies -
- * The force is measure in kg/cm² -

Result -

- Paracetamol tablets was prepared & evaluated. -

Object:-

~~ABSTRACT~~ evolution of glass container as per I.P.-

Reference-

"DR. shekh Aijaz, DR. kailash" industrial pharmacy-I
practical book published by Pee Vee.- Edition : 2018 &
page number : 95 to 99 -

Requirements-(A) Apparatus-

• Beakers, measuring cylinder,
measuring cylinder, & autoclave etc-

(B) Chemicals-

• Acetone & sulphuric acid -

Theory-

- Glass is employed as the container material of choice for most of pharmaceutical products.-
- Glass containers are evaluated by USP tests powdered glass test & water attack test-

* Types of glass-

- Type-I : Boro-silicate glass-
- Type-II : Sodalime treated glass-
- Type-III : Sodalime glass-
- Type-IV : NP general purpose sodalime glass not for parenterals-

Procedure-

* Powdered glass test-

• Step-1: Preparation of glass specimen-

Few containers rinsed thoroughly with purified water & dried with stream of clean air, grind the containers in a mortar to a fine powder & pass through sieve no. #20 & #50.-

• Step-2: Washing the specimen-

10gm of the above specimen is taken into 250ml conical flask & wash it with 30ml acetone.-

• Repeat the washing, decant the acetone & dried after which is used within 48hrs.-

Powder glass test -

Types of glass	Amt. of titrant consumed (ml)
• Type - I	0.1 ml
• Type - II	8.5 ml
• Type - V	15 ml

Date...15/10/22...

S. No...11...

Page No...33...

• Step-3 -

- 10 gm of sample is added with 50ml of high purity water in 250ml flask.
- Place it in an autoclave at 121°C for 30min.
- Cool it under running water. Decant the solution into another flask, wash again with 15ml distilled water & again decant.
- Titrate immediately with 0.02N sulphuric acid using methyl red as an indicator & record the volume.

Result :-

- Glass containers were passed the test as per IP.